**FR Informations et instructions concernant le produit / EN Product Information and Instructions / IT Informazioni sul prodotto ed istruzioni per l'uso / DE Produktinformationen und Gebrauchsanweisung / ES Instrucciones e información sobre el producto / NL Productinformatie en instructies / DA Produktinformation og -vejledning / FI Tuotetta koskevat tiedot ja käyttöohjeet / SV Produktinformation och instruktioner / PT Informações do produto e instruções / EL Πληροφορίες προϊόντος και οδηγίες / PL Informacje o produkcie oraz instrukcje / CS Informace o výrobku a návod k použití / HU A termék ismertetése és használat / TR Ürün Bilgileri ve Talimatları / NO Produktinformasjon og instruksjoner**

**Instruments Monopolaires pour Electrochirurgie – Tubes métalliques avec gaine amovible (Utilisation laparoscopique) / Monopolar Electrosurgical Instruments – Metallic tubes with removable sheath (Laparoscopic Use) / Strumenti elettrochirurgici monopolari -Tubi metallici con guaina rimovibile (per uso laparoscopico) / Monopolare elektrochirurgische Instrumente - Metallrohre mit abnehmbarer Schleuse (für Laparoskopie) / Instrumental electroquirúrgico monopolar - Tubos metálicos con vaina extraíble (para cirugía laparoscópica) / Monopolaire elektrochirurgische instrumenten - Metalen buizen met verwijderbare sheath (voor laparoscopisch gebruik) / Monopolære elektrokirurgiske instrumenter - Metalrør med aftageligt hylster (laparoskopisk brug) / Monopolaariset sähkökirurgiset instrumentit - Metalliputket, joissa on irrotettava holkki (laparoskooppinen käyttö) / Monopolära elektrokirurgiska instrument - Metallrör med avtagbar hylsa (laparoskopisk användning) / Instrumentos electroquirúrgicos unipolares - Tubos metálicos com bainha removível (utilização laparoscópica) / Μονοπολικά ηλεκτροχειρουργικά εργαλεία - Μεταλλικοί σωλήνες με αφαιρούμενο θηκάρι (Λαπαροσκοπική χρήση) / Jednobiegunowe przyrządy elektrochirurgiczne - Metalowe rurki ze zdejmowalną osłoną (zastosowania laparoskopowe) / Monopolární elektrochirurgické nástroje - Kovové trubice se snímatelným sheathem (laparoskopické použití) / Monopoláris elektrosebészeti eszközök - Fémcsővek levehető hüvellyel (laparoszkópos használatra) / Monopolar Elektro Cerrahi Aletleri - Çıkarılabilir kilifli metalik tüpler (Laparoskopik Kullanım) / Monopolare elektrokirurgiske instrumenter - Metallrør med avtakbar hylse (laparoskopisk bruk)**

**Manufacturer:**  
Integra MicroFrance SAS  
Le Pavillon • 03160 Saint Aubin Le Monial • France  
integralife.com

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## EN English / Product Information and Instructions

### DESCRIPTION OF THE DEVICES

Monopolar electrosurgical instruments are available in different configurations to be used in laparoscopic and open surgery as indicated by surgical techniques. Electrical insulation is applied to selected portions of the instrument and a connector is provided to attach the instrument to a variety of currently available monopolar electrosurgical generators. These devices are intended to be used by qualified physicians and operating room personnel familiar with high frequency current and who have been properly and correctly trained in the surgical techniques that must be used. Physicians should seek relevant pre-clinical education before attempting new surgical procedures, especially those in insufflated cavities.

## **MATERIALS**

The instruments are fabricated of one of the following biocompatible materials: glass, ceramic, titanium nitride, stainless steel, tungsten, thermoset polymers (including silicone), thermoplastic polymers, sterling silver or chromed brass.

## **INTENDED USE**

Monopolar electrosurgical instruments are indicated for laparoscopic and non-laparoscopic procedures where the physician determines that high-frequency electrical current is appropriate to remove tissue and/or control bleeding.

## **INDICATIONS FOR USE**

The electrosurgical instruments are intended to remove tissue and/or control bleeding.

## **CONTRAINDICATIONS**

None known.

## **SAFETY**

- Before each use, check the outer surface of the parts introduced into the patient to verify that there is no rough surface, no sharp edge or undesirable protrusion which may cause a safety hazard.

## **WARNINGS**

- The maximal assigned output voltage that can be applied to a monopolar high-frequency electrosurgical accessory is 1355 Vrms [1958 Vpeak] except for special indications etched on the instrument and/or on a specific insert. Such peak voltage can be achieved with certain modes of functioning of some electrosurgical generators. Please check with the manufacturer of the electrosurgical generator and/or refer to the technical documentation supplied by the manufacturer of the electrosurgical generator.
- Do not use these instruments for purposes other than those for which they are intended.
- Do not bend, disassemble or force an instrument; breakage or failure of the instrument could occur, resulting in possible harm to the patient or user.
- Use appropriate personal protection equipment when handling contaminated instruments.
- For devices that incorporate suction, do not extend suction tube while monopolar energy is applied, in order to avoid accidental burns.
- Ensure that the tissue intended to be coagulated is grasped in the jaws before applying monopolar energy, to avoid harming tissue.
- Do not energize the electrosurgical instrument in the presence of flammable anesthetics. Fires can occur if this type of instrument is used in an environment rich in oxidizing gases such as oxygen or nitrous oxide (N<sub>2</sub>O), including intubation tubes (burns of the respiratory pathways) and the concentrations of oxygen collecting under surgical drapes. Review security procedures before using this instrument.
- Non-flammable agents should be used for cleaning and disinfection wherever possible. Flammable agents used for cleaning or disinfecting or as solvents of adhesives should be allowed to evaporate before the application of HF surgery. There is a risk of pooling of flammable solutions under the patient or in body depressions such as the umbilicus or in body cavities such as the vagina. Any fluid pooled in these areas should be mopped up before using the instrument. Attention should be called to the danger of ignition of endogenous gases. Some materials, for example cotton or gauze, when saturated with oxygen, may be ignited by sparks produced in normal use.
- Do not place the instrument on the patient when it is not being used. Place the instrument in an insulated support or on a clean, dry surface, very visible and non-conductive, so as to avoid accidental electrical injuries.
- Do not try to modify the instrument. Do not try to repair the electrical insulation.
- Do not "buzz" the instrument during the surgical procedure to reduce the risk of burning the patient or the physician.
- To avoid lesions at other sites, only activate the generator when the surgical field and the distal extremity of the forceps are sufficiently visible.
- Use extreme care when coagulating long, thin tissue, such as adhesions. These tissues can carry electrical current to an unseen or remote location. In this case, the use of bipolar techniques can be desirable so as to avoid tissue damage.
- Trocar canulas, if used with the device, should be inspected for sharp edges that may damage the electrical insulation of the instrument. When using a trocar canula, the working end of the device must pass entirely through the trocar canula for the device to function properly and effectively.
- The pathways of the current through conductive elements like metal instruments and endoscopes can cause local burns to the patient, the physician or another member of the care team. Contacting conductive elements with the active cautery area may cause undesired tissue heating and burns.
- The temperature of the instrument in the active cautery area can remain high enough after use that it burns the patient, the physician or another person, even when the electrical current is turned off.
- Alternate site injuries may occur at an electromyography (EMG) needle electrodes monitoring site.
- Electrosurgical devices may cause electromagnetic interference with other operating room devices.
- To avoid alternate site burns, ensure hinges, handles and proximal un-insulated portions of the electrosurgical instruments do not inadvertently contact patient.
- For patients with cardiac pacemaker or other active implants, a possible hazard exists because interference with the action of the pacemaker may occur or the pacemaker may be damaged. In case of doubt, advice or approval from the manufacturer of the active device should be obtained.
- Extreme care should be taken when handling instruments with a dielectric coating. Damage to the dielectric coating may result in patient / user injury.
- Begin procedure at the lowest possible electrosurgical power setting to reduce the risk of patient burns at high voltages.
- To ensure proper functioning of dismantable devices and devices with accessories, check the assembly and the functionality of all elements of the device before use.
- Connect the device only to the monopolar receptacle of the generator. Improper connection may result in inadvertent activation or other potentially hazardous conditions.
- The device requires the use of monopolar electrosurgical generators that monitor the current of the return electrode.
- Apply the patient return electrode according to the recommendations of the generator manufacturer.
- The entire surface of the neutral electrode should be securely connected to the patient's body and as close as possible to the surgical field.
- The patient should not be in contact with grounded metal parts or parts having an appreciable capacity with respect to the ground (for example operating table, supports, etc.). Antistatic wrapping is recommended in this case.
- Skin to skin contact (for example between the patient's arms and body) must be avoided, for example by separation with dry gauze.
- Discard instrument after suspected Creutzfeldt-Jakob Disease (CJD) exposure; Integra instruments have not been validated to withstand the chemical and thermal exposure recommended to eradicate prions.

## **PRECAUTIONS**

- Do not use abrasive cleaning products such as hard brushes, etc.
- Make a visual inspection of the instrument and the cable to ensure that the electrical insulation is in good condition.
- Extreme care must be taken when handling and cleaning instruments with sharp points or edges.
- The use of trocars made entirely of plastic or metal is recommended with electrosurgical devices in order to avoid thermal injuries in the surgical access zones.
- We recommend you to use MICROFRANCE® trocars.
- Delicate instruments should be handled with extreme care to prevent damage.
- Limit coagulation / cutting with scissors as much as possible to preserve the blades.

## **ADVERSE EVENTS**

- Localized burns to the patient or physician may result from electrical current carried through conductive objects (such as trocar cannulas). Electrical current may be generated in conductive objects by direct contact with the active electrode, or by the active accessory (electrode or cable) being in close proximity to the conductive object.

#### INSTRUCTION FOR USE

- The instrument is intended to be used with electrosurgical generators and cables compliant with standards IEC60601-2-2 and/or AAMI HF-18 using a neutral electrode (REM).
- Reusable instruments must be cleaned and sterilized before the first use and before every new use.
- It is the responsibility of the surgical team to select the appropriate instruments for each use.
- Ensure the electrosurgical generator is in off or standby mode, and connect the device cable to the monopolar receptacle of the electrosurgical generator only, ensuring the connector pins are fully seated in the generator receptacle.
- Position electrosurgical cables to avoid contact with the patient, user or other cables.
- Begin procedures at the lowest possible electrosurgical power setting. Slowly increase power until the desired coagulation effect is achieved. This practice reduces the possibility of the following problems: a) capacitive coupling b) damage to insulation c) increased risk of patient burns at high voltages.
- If the output power appears weak at normal functional settings, check connections, cables, the generator and the instrument before selecting higher output power.
- Carefully insert and remove the instruments from the trocars. Close the jaws of the instruments when passing through trocars.
- For instruments that incorporate suction, a stylet is recommended clearing and/or removing clogged tissue.
- Avoid prolonged activation. Short, multiple activations prevent surrounding tissue from heating up.

#### REPROCESSING INSTRUCTIONS for Reusable Instruments (PER ISO17664)

When applying dry times to Integra instrument trays and their accessories, dry times outside the standard healthcare prevacuum parameters may be required. The current recommended dry times may range from a standard 15 minutes to an extended time of 40 minutes. This is especially important for polymer-based (plastic) instrument trays. The user should have adequate means to verify drying time, such as visual inspection.

<b>Warnings/ Precautions</b>	To prevent stains, use distilled or demineralized water, and uses a neutral detergent to reprocess these instruments. Do not cold soak in glutaraldehyde, chlorine, or ammonium solutions, or dry heat sterilize, as damage to the instrument finish may occur.																						
<b>INSTRUCTIONS</b>																							
<b>1) Bringing into service of a new instrument</b>	This product is provided non-sterile and must be cleaned and sterilized before the first use and any reuse. The protection / silicone caps used during shipment must be removed before cleaning and sterilization procedures.																						
<b>2) Preparation at the point of use</b>	It is recommended that instruments be reprocessed as soon as is practical following use. Rinse instrument immediately after use in distilled or demineralized water. Keep instrument moist (for example, cover with a wet drape or saturate with tap water – do not use saline) after use and prior to decontamination, cleaning and sterilization to ensure adequate cleaning. Instruments with cleaning ports and/or suction tubes shall be flushed with neutral pH enzymatic detergent, followed by a tap water rinse.																						
<b>3) Preparation before cleaning</b>	Fully disassemble modular designed instrument for effective cleaning. Remove any cap covering the cleaning port, if applicable. For instruments with lumens and instruments equipped with cleaning ports, inject cleaning solution through the instrument with an irrigation syringe. When cleaning suction tubes with holes to control suction, place gloved finger over the control hole to flush through the tube. Open any articulated instruments before positioning in the soaking solution. Place instruments so they do not touch each other. Note: Flushing in running water is essential between decontamination and cleaning to prevent any risk of reaction between the two solutions.																						
<b>4) Cleaning, disinfection and drying</b>	<p><b>Cleaning: Automated:</b> All instruments shall be thoroughly rinsed with tap water prior to placement in the automatic washer. Remove instruments and equipment from any sterilization trays before placing into washer baskets. Orient devices following recommendations of washer/disinfector manufacturers. Use alkaline or neutral pH detergent recommended by washer/disinfector or detergent manufacturers. These products have been validated for effective cleaning using an automatic washer/disinfector cycle consisting of the following:</p> <table border="1" data-bbox="363 1305 1509 1509"> <thead> <tr> <th>Phase</th> <th>Recirculation Time</th> <th>Water temperature</th> <th>Detergent type</th> </tr> </thead> <tbody> <tr> <td>Pre-wash 1</td> <td>3 minutes</td> <td>cold tap water</td> <td>N/A</td> </tr> <tr> <td>Enzyme Wash</td> <td>5 minutes</td> <td>hot tap water</td> <td>neutral pH enzymatic detergent</td> </tr> <tr> <td>Wash 1</td> <td>5 minutes</td> <td>66°C set point</td> <td>neutral pH detergent</td> </tr> <tr> <td>Rinse 1</td> <td>3 minutes</td> <td>hot tap water</td> <td>N/A</td> </tr> </tbody> </table> <p>Thoroughly examine instruments for any residual soil. A rigid MIS (minimally invasive surgery) rack should be used to clean rigid MIS instruments and instruments with lumens to maximize efficiency of cleaning. Suction ports and flush ports should be connected to the MIS Rack to allow a constant flow of solution inside hard-to-reach areas of the instruments.</p> <p><b>Cleaning: Manual:</b> Soak in lukewarm (less than 43°C), mild (pH 7.0 – 8.5), enzymatic detergent and deionized water for a minimum of 2 minutes. For instruments with lumens and instruments equipped with cleaning ports, inject cleaning solution through the instrument with an irrigation syringe and allow to soak for a minimum of 2 minutes. Use a soft instrument brush to scrub instruments while submerged in cleaning solution to remove organic matter. Rinse with deionized water, then clean ultrasonically in a lukewarm (less than 43°C), mild (pH 7.0 – 8.5) detergent and deionized water for 10 minutes. Rinse thoroughly with deionized water, utilizing a syringe to thoroughly rinse cleaning solution from lumens and cleaning ports. Clean in this manner until no visible soil remains on the instrument. Dry with compressed air or wipe dry with a lint-free cloth. Examine instruments for any staining or deterioration, remove from use as appropriate. Note: When using an ultrasonic cleaner or a spray washing machine, follow the manufacturer's recommendations, particularly with regard to articulated instruments and positioning of instruments. Note: Additional cleaning methods may be warranted, including presoaking in 3% hydrogen peroxide.</p>			Phase	Recirculation Time	Water temperature	Detergent type	Pre-wash 1	3 minutes	cold tap water	N/A	Enzyme Wash	5 minutes	hot tap water	neutral pH enzymatic detergent	Wash 1	5 minutes	66°C set point	neutral pH detergent	Rinse 1	3 minutes	hot tap water	N/A
Phase	Recirculation Time	Water temperature	Detergent type																				
Pre-wash 1	3 minutes	cold tap water	N/A																				
Enzyme Wash	5 minutes	hot tap water	neutral pH enzymatic detergent																				
Wash 1	5 minutes	66°C set point	neutral pH detergent																				
Rinse 1	3 minutes	hot tap water	N/A																				
<b>5) Maintenance, inspection and testing</b>	Check the cleanliness and operation of the instrument. Clean again if debris is present and remove from use any damaged instrument. Inspect components for any damage. If damage is observed, do not use the instrument until it is repaired. Following cleaning, lightly lubricate instruments with movable parts. Use a lubricant intended for sterilizable instruments such as a water-soluble instrument milk. Do not use silicone spray.																						
<b>6) Packaging</b>	A standard sterilization wrap may be used. In the US, an FDA cleared surgical wrap is required. In sets: Instruments may be loaded into dedicated instrument trays or general purpose sterilization trays. Ensure that cutting edges are protected. Wrap trays using appropriate method. <b><u>The complete moistening of lumen instruments, prior to sterilization, is required for an optimal efficiency of the steam sterilization cycle.</u></b>																						

<b>7) Sterilization (Temperatures are minimum required; times are minimum required)</b>	The sterilization parameters given below should be used for devices that are fully disassembled when disassembly is possible. During sterilization, do not mix instruments of different metals, to avoid electrolyte oxidation. Close instruments with catches and racks on the first notch. Arrange the instruments in sterilization containers with perforations on the top and bottom and on supports such as those used in microsurgery. Follow the appropriate cycle listed in the table below. All steam cycles have been validated in the wrapped configuration and can be sterilized wrapped or unwrapped. These devices have only been validated for steam sterilization methods.					
		Instruments Only		Instruments Only or Instruments in Instrument Trays		
	Cycle	Gravity	Gravity	Pre-vac	Pre-vac (FR/WHO)	Pre-vac (U.K.)
	Temperature	121 °C	132 °C	132 °C	134 °C	134 °C
	Time	30 min	10 min	4 min	18 min	3 min
<b>8) Storage</b>	Store instruments in a clean, dry area.					
<b>9) Disposal</b>	Instruments must be disposed in accordance with applicable local regulations.					

The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. They are not applicable to single use devices or single use accessories, which must be destroyed after use in accordance with applicable local regulations. It remains the responsibility of the processor to ensure that reprocessing is performed using validated equipment to achieve the desired result. This normally requires validation and routine monitoring of the process. Some devices have specific assembly instructions. In this case, refer to the assembly insert provided with the device for additional instructions. All validations performed per AAMI TIR12, (Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers).

## RETURNS

Contact Integra Customer Service to obtain a Return Goods Authorization number (RGA #) prior to shipping the product back to Integra. The RGA # should be prominently displayed on the outside of the return packaging and included on all paperwork enclosed with the return. All product returned should be decontaminated and cleaned, and all products should be safely packed in protective wrapping. Please have the original invoice number or purchase order number available to assist in verifying warranty information.

## LIMITED WARRANTY

A. This LIMITED WARRANTY provides assurance for the customer who purchases a Integra Product (hereinafter the "Product") that should the Product fail to function to Integra's published specifications during the term of this LIMITED WARRANTY (one year from the date of shipment for new Product), Integra will either replace, repair, or issue a credit (adjusted to reflect the age of the Product) for the Product or any portion thereof. This LIMITED WARRANTY is extended only to the buyer purchasing the Product directly from Integra or from its affiliate or its authorized distributor or representative.

B. To qualify for this LIMITED WARRANTY, the following conditions must be met:

- (1) The Product must be used on or before its "Use By" or "Use Before" date, if applicable.
- (2) The Product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.
- (3) Integra must be notified in writing within thirty (30) days following discovery of a defect.
- (4) The Product must be returned to Integra within thirty (30) days of Integra receiving notice as provided for in (3) above.
- (5) Upon examination of the Product by Integra, Integra shall have determined that: (i) the Product was not repaired or altered by anyone other than Integra or its authorized representative, (ii) the Product was not operated under conditions other than normal use, and (iii) the prescribed periodic maintenance and services have been performed on the Product.

C. This LIMITED WARRANTY is limited to its express terms. THIS LIMITED WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED WHETHER STATUTORY OR OTHERWISE, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. In no event shall Integra be liable for any consequential, incidental, prospective or other similar damage resulting from a defect, failure, or malfunction of the Product, whether a claim for such damage is based upon the warranty, contract, negligence or otherwise.

D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. Users may benefit from statutory warranty rights under legislation governing the sale of consumer goods. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.

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